

stryker



CORE 2 console

**Powering**  
your instruments

CORE 2 console

# Driving benefits



## Versatility

- For use in a variety of surgical procedures, including orthopaedic, dental, ENT, neuro, spine, and endoscopic applications
- Powers Stryker's small bone saws and drills, NS high speed drills, ENT shavers, joint shavers and bone mill

## Unparalleled customization

- 100+ user profiles that can be stored and transferred
- A Adjustable acceleration and braking**
- Adjustable irrigation flow rate
  - Largest breadth of compatible product portfolios
  - Copy/edit allows creation of multiple profiles without starting from scratch
  - Customizable features instantly appear based on the device connected and include RPMs, oscillation rate, direction, fixed/variable mode and many more
  - Transferable to other CORE 2 consoles via USB

## Safety

- Built-in design features make safe operation as effortless as possible
- B Constant, color-lit visibility of matched drills and foot pedals**
- 13-year proven performance
  - Hand piece (C) and foot pedal (D) ports now flush with surface, not recessed

## Workflow efficiency

- Powers a spectrum of devices, reducing need for multiple power sources
- C Three hand piece ports; runs two hand pieces simultaneously**
- B Colored illumination rings (5) display matched drills, saws and foot pedals for free and confident workflow**
- D Two foot pedal ports; ability to assign two pedals to a single hand piece**
- E Integrated irrigation; one-step insertion and removal**

## Simplicity

- Smart programming auto-activates screens, options and information based on connected devices
  - Large icons and 7" color touch screen with 170° wide viewing angle
  - **Locations of hand piece (C) and foot pedal (D) have been flipped to provide more room around irrigation cassette port**
  - Numeric and written error descriptions with troubleshooting tips
  - Flat, intuitive navigation with shortcuts for added speed
  - Stryker sales rep contact info can be stored internally
- F Distinct audio tones denote different actions/alerts**

# Core solution

CORE stands for Consolidated Operating Room Equipment, but its figurative meaning rings true as well. Our console is core to what you do. Paired with our market leading saw and drill portfolios, this console works across multiple surgical specialties. Here's a snapshot of how the CORE 2 console plays an integral role in optimizing your OR investments and efficiency.

## Leading drills and saws, leading console

Our console is designed with features and software that help you gain optimal benefits from our saws and drills, including those shown here.

TPX micro drill



TPX oscillating saw



TPX sagittal saw



TPX reciprocating saw



TPX universal driver



TPX wire driver



## Order or learn more

Here are just some of the CORE 2 console components available.

For a complete listing or more information please call your Orthopaedic Instruments sales representative, **800 253 3210** or visit [stryker.com/us/en/orthopaedic-instruments.html](http://stryker.com/us/en/orthopaedic-instruments.html)

**CORE 2 console**  
5400-052-000US

**Single handpiece irrigation cassette (disposable)**  
5400-050-001

**Dual handpiece irrigation cassette (disposable)**  
5400-050-002

**TPX handpiece cord**  
7400-004-000

**CORE accessory cart**  
5400-410-000

**Bi-directional footswitch**  
5100-008-000

**Uni-directional footswitch**  
5100-007-000

### TPX handpieces

**TPX micro drill**  
7400-015-000

**TPX oscillating saw**  
7400-031-000

**TPX XL oscillating saw**  
7400-131-000

**TPX sagittal saw**  
7400-034-000

**TPX reciprocating saw**  
7400-037-000

**TPX universal driver**  
7400-099-000

**TPX wire driver**  
7400-062-000

**TPX handswitch**  
7400-009-000

## Your system. Your way.

With foot pedals, irrigation cassettes, a mobile cart and many other accessories, you can create the system you want and make its management easy.

### First-rate foot pedals

Get your choice of two foot pedals, each with multiple buttons to set desired functions such as irrigation on/off, forward/reverse direction, change drill port assignment, change RPMs and more.



**Bi-directional footswitch**

Dual pedal with three feature buttons

**Uni-directional footswitch**

Single pedal with two feature buttons



### No fuss irrigation

Our front facing, integrated irrigation with one-step insertion and removal makes set up simple, without any extra hardware components to affix or tubes to thread. Users also love its integrated I.V. pole bracket, auto prime and adjustable flow rate.

### Convenient cart

With eight drawers, two baskets and a top surface area, this cart holds your CORE 2 console, foot pedals and a full inventory of cutting accessories. Organize as you wish and apply custom drawer labels to make management easy.



## Orthopaedic Instruments

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate a Stryker product. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: CORE, Stryker, TPX. All other trademarks are trademarks of their respective owners or holders.

D0000040324 AB  
Copyright © 2021 Stryker  
Printed in USA

Stryker  
1941 Stryker Way  
Portage, MI 49002 USA  
t: 269 323 7700 f: 800 999 3811  
toll free: 800 253 3210

[stryker.com/us/en/orthopaedic-instruments.html](http://stryker.com/us/en/orthopaedic-instruments.html)