

stryker

Innovation that **gets to the heart**



System 9
SABO
sagittal
saw

System 9
Sternum
saw

Unmatched performance, whenever you're ready

We innovate constantly to enhance your experience without changing the core of what gives our saws the most important power of all—the bond between a surgeon and the tool in their hand.

- ✔ Well-balanced and designed to fit in your palm while reducing the number of pressure points while cutting
- ✔ Over time, exposure to pressurized water and harsh chemicals can wear down the toughest power tools—until now. Introducing our most robust anodized coating paired with IPX 9 rating for lasting performance

SABO sagittal saw

Our sagittal saw stands up to every case a surgeon takes on

- ✔ Blade mount is easy to reposition for a more accurate cutting experience
- ✔ Intentional design allows you to cut closer to the surgical site*
- ✔ Small lightweight handpiece for improved control

Sternum saw

Our sternum saw is held to a different standard—yours

- ✔ Sternum blade guards protect soft tissue from coming into contact with blade
- ✔ Increased versatility with two cut speed options: 11,000 or 14,000 cuts per minute
- ✔ A wide variety of sternotomy blades available



* in comparison to the cordless driver with sagittal saw attachment

Batteries

Efficiency through tactile and visual feedback

- ✔ The new handpiece-to-battery interface provides tactile and visual feedback, enabling efficient assembly and confident connection
- ✔ The battery's 45-day shelf life reduces risk of compromised performance
- ✔ System 9's battery indicators are designed to be brighter and repositioned, creating greater visibility for surgeons and the surgical team
- ✔ The battery housing features new contacts so that your battery is more dependable and durable



Wireless charger

Our new wireless charger is the first of its kind—an advanced piece of technology for a smarter, more streamlined OR experience

- ✔ Six bays to place your batteries—and your trust
- ✔ Most versatile battery charging capabilities to date, allowing you to charge your batteries in any environment that accommodates your workflow
- ✔ Ability to confirm battery status pre-procedure
- ✔ Simplified user interface
- ✔ Ability to charge battery post sterilization while maintaining sterility, including ability to charge through blue wrap



Get started

For questions or to place your order, please contact your sales rep or contact Customer Service.

Power tools and accessories

Product no.	Description
4608-000-000	SABO sagittal saw
9207-000-000	Sternum saw
7207-003-000	Sternum blade guard
7207-002-000	Sternum blade guard long

To order

Customer Service | P 800 253 3210 | F 877 790 7772 | inst.stryker.cs@stryker.com

Orthopaedic Instruments

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: Flex Financial, ProCare, SABO, SEM, Smart Equipment Management, Stryker and System 9. All other trademarks are trademarks of their respective owners or holders.

D0000257996 Rev. AA | Copyright © 2023 Stryker

[stryker.com](https://www.stryker.com)